



**Artwork and Signature File for:**

**MAN-00850, "IHE INTEGRATION STATEMENT, SELENIA DIMENSIONS"**

Artwork consists of:

- One (1) 8 ½ inch x 11 inch sheet(s) attached.

REV AUTHORED BY	DATE				
J. SHARP	12/16/09				
REV DRAFTED BY	DATE	TITLE		DOCUMENT NUMBER	REV
J. SHARP	12/22/09				
PROPRIETARY INFORMATION: The content of this document is the exclusive property of Hologic and may not, without prior written permission of Hologic, be reproduced, copied or used for any purpose whatsoever.		AW, IHE INTEGRATION STATEMENT, SELENIA DIMENSIONS		AW-02185	002
		ARTWORK		SIZE A	SHEET 1 OF 1

**IHE Integration Statement for Selenia Dimensions**

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product</b>	<b>Version</b>	<b>Date</b>
Hologic, Inc.	Selenia Dimensions Acquisition Workstation (AWS)	1.0.x	December 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow	Acquisition Modality	Assisted Acquisition Protocol Setting	
		PPS Exception Management Option	
		Broad Worklist Query	
		Patient Based Worklist Query	
Patient Information Reconciliation	Acquisition Modality	None	
Mammography Image Profile	Acquisition Modality	None	
Mammography Image Profile	Print Composer	None	
Access to Radiology Information	Image Display	None	
Portable Data for Imaging	Portable Media Creator	None	
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	Not applicable		
<b>DICOM</b>	<a href="http://www.hologic.com/product-support/digital-mammography/dimensions/">www.hologic.com/product-support/digital-mammography/dimensions/</a>		
<b>Links for general information on IHE</b>			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	