

Artwork and Signature File for:

MAN-01288 “IHE INTEGRATION STATEMENT, SEL 3.4”

Artwork consists of:

- One (1) 8 ½ inch x 11 inch sheet attached.

REV AUTHORED BY	DATE	HOLOGIC™	SIGNATURES ON FILE	
A. TAMBASCIO	03/27/09			
REV DRAFTED BY	DATE	TITLE IHE INTEGRATION STATEMENT FOR SELENIA V3.4	DOCUMENT NUMBER	REV
J. SHARP	04/01/09		AW-03168	001
PROPRIETARY INFORMATION: The content of this document is the exclusive property of Hologic and may not, without prior written permission of Hologic, be reproduced, copied or used for any purpose whatsoever.			ARTWORK	SIZE A

IHE Integration Statement for Selenia v3.4.x

IHE Integration Statement			
Vendor	Product	Version	Date
Hologic	Selenia Full Field Digital Mammography Acquisition Workstation (AWS)	3.4.X	March 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Assisted Acquisition Protocol Setting	
		PPS Exception Management Option	
		Broad Worklist Query	
		Patient Based Worklist Query	
Patient Information Reconciliation	Acquisition Modality	None	
Mammography Image Profile	Acquisition Modality	None	
Access to Radiology Information	Image Display	None	
Links to Standards Conformance Statements for the Implementation			
HL7	Not applicable		
DICOM	http://hologic.com/cc/dicom.htm		
Links for general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	